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English

1. Product Function

Experience gained during the last 30 years in cruciate ligament surgery has taught that the natural cruciate ligament best reflects the physiological kinematics and dynamic stabilization of the knee. Following this principle, the technique of Dynamic Intraligamentary Stabilization (DIS) after injury of the anterior cruciate ligament was developed.

Using the Ligamys implant, the knee joint is stabilized in a controlled manner by a dynamic spring system. The local ligament receives the immobilization necessary for stable scarring while simultaneously being exposed to a dynamic stimulus that ensures better structuring of the healing ligament. The spring system developed in the context of these considerations, which is used for precise reconstruction of the dynamics of the ligament, provides both functions: Stability and dynamics.

The dynamic system (DIS-Ligamys), consisting of a sleeve with a spring system and an UHMWPE yarn with a «button», makes it possible for the knee joint to be artificially kept in the posterior drawer during the first phase of healing, enabling and supporting scarring of the ligament stumps.

Only natural scarring healing of the tissue ensures functional restoration.

2. General Warnings and Safety Notes

Before using the Ligamys implant manufactured by Mathys Ltd Bettlach, users must familiarise themselves with the handling of the instruments and the product-specific surgical technique as well as with the warnings, product-specific risks, safety notes and recommendations listed in the package insert. Compliance with the training concept and surgical technique specified by Mathys Ltd Bettlach is mandatory. For information, please contact the competent national representative.

The implant may be used only by a physician trained by Mathys Ltd Bettlach.

a) Notes on Storage and Application

Ligamys implants are designed to be used only once. All implants are individually packaged under sterile conditions in protective packaging. The implant must be carefully stored in its unopened protective packaging. This protects the implant from external influences. The storage conditions stated on the packaging must be observed. Store the implant in a dry place at room temperature. The expiration date indicated on the packaging must be observed. After the expiration date, the implant must not be used anymore. Before unpacking the implant, inspect the protective pack-

aging for damage. The protective packaging must not be opened until immediately before use. Subsequently, the implant is to be implanted immediately after removal from the packaging. When removing the implant from the sterile package, observe the relevant regulations concerning aseptic techniques.

Prior to implantation, all components of the implant must be checked for potential damage. Only undamaged implants may be used. Modular components are integral parts of a system and may therefore be combined only with appropriate system components. Re-use of an implant previously inserted is prohibited. Implants and system components may be neither mechanically processed nor modified.

b) Notes on Sterility

The package is marked with the word «STERILE» and consists of the protective packaging and a sterile packaging. Sterility is ensured only while the sterile packaging is undamaged. Implants marked as «STERILE» have been sterilized using ethylene oxide (EtO) or β/γ irradiation. In case of EtO sterilization, there is an appropriate marking on the label.

c) Notes on Information Supplied

Before any surgery, make sure the implant matches the information on the package and in the package inserts (article number, identification number, product name). For further information (product name, item

number, batch number, etc.) please see the package, the included adhesive labels and the package inserts. For detailed information on the material of the components used, please refer to the relevant pages of the product catalogue or the corresponding product information. The adhesive labels included must be stuck into the patient record and into the surgical report. Thus, traceability of the implant is always ensured.

d) Pre-, Intra- and Post-operative Warnings and Safety Notes

Before using the products, users must familiarize themselves with the general issue of rupture of the anterior cruciate ligament (ACL) and the product-specific surgical technique. The recommended surgical technique is to be complied with. Patient selection greatly depends on the patient's age, the age of the primary anterior cruciate ligament rupture, general health, condition of the existing bone material, previous surgery and anticipated additional surgery. Patients should be informed about the following with regard to product, product-specific surgical procedure and product-specific information on rehabilitation:

- The risks of implantation as well as the possible impacts of the conditions named under items 3 to 5 and the possible complications mentioned in item 6
- The necessary post-operative treatment (particularly concerning possible activities and weight loads)
- The need for any unusual changes at the operated site to

be reported to the treating physician

e) Additional Warnings and Safety Notes

Responsibility for selection and use of the implant rests exclusively with the user. Mathys Ltd Bettlach or your local representative will take the implant back only if the protective packaging is unopened and fully intact. Mathys Ltd Bettlach does not guarantee sterility if the protective packaging is damaged or the sterile package has been opened improperly. The packaging material of opened implant packages must be disposed of or returned to the manufacturer for disposal. Reprocessing of the implant is not permitted and releases the company from all liability.

3. Indications

- Fresh primary rupture of the anterior cruciate ligament, not older than 21 days on the day of surgery.

4. Contraindications

- Acute or chronic infections, local or systemic (or corresponding history)
- Severe muscle, nerve or vascular disease that may endanger the affected limb
- Hypersensitivity to the materials used (e.g., cobalt, chromium, nickel, etc.)
- Insufficient bone substance, poor bone quality or insufficiently mature bone, which might compromise stable anchorage of the implant
- Any circumstances that may prevent the patient from appropriately limiting his/her activities or following

medical instructions during the healing phase

- Patients in whom conservative therapy is promising

5. Risk Factors

- This product may be unsuitable for patients with inadequate or immature osteogenesis. In this patient population, the placement of the implant must not affect the epiphyseal plate
- Lack of compliance
- Bone disease
- Severe deformities
- Infection or generally weakened immune system
- Drug or alcohol abuse
- Hypersensitivity to the material to be implanted

6. Possible Complications

- General complications associated with surgery in general, drugs, instruments additionally used, blood, etc.
- Specific implant-associated complications that may arise from the factors mentioned under item 5 (risk factors) (e.g. shifting of the implant components, fatigue failure, early or late loosening, abrasion or wear, etc.)
- Pain
- Tissue reactions caused by allergy or foreign body reactions to abraded particles
- Loosening and migration of the implant
- Haematomata and delayed wound healing
- Temporary or permanent neuropathy
- Venous thrombosis or pulmonary embolism
- Postoperative infection with possible sepsis

7. Mathys Ligamys implants in the context of magnetic resonance (MR)

The safety rating «MR Conditional» means that in

a specific MR environment with defined application conditions a product demonstrably poses no known hazards.

Notes on MRI Safety

Non-clinical tests have shown Ligamys implants to be MR Conditional. A patient with such an implant can be safely examined in an MR system that satisfies the following conditions:

- Static magnetic field strength of 1.5 or 3.0 Tesla, respectively
- Maximum spatial gradient field of 3000 Gauss/cm (30 T/m)
- Maximum specific absorption rate specified for the MR system, whole body average (SAR), of 4 W/kg for 15 minutes of scanning

Note: The defined SAR must be strictly adhered to if the nearest part of the implant is closer than 30cm to the isocentre of the MR scanner. If the implant is in its entirety more than 30cm away from the isocentre of the MRI device, there is lower RF-induced heating of the tissue in the vicinity of the implants, and MRI scans in normal operating mode are permissible for all implants.

- Normal mode of the MR system. The effects of MRI procedures with MR systems and conditions that exceed these values have not been determined.
- The impact of local RF transmission coils has not been studied, and these are not recommended in the area of the implant.

Attention: Combination with other products has not been evaluated for safety and compatibility in the MR environment. It has not been tested for RF heating in the MR environment. It is not known whether combination with other products is safe in the MR environment. Examination of a patient with such a product combination may result in patient injury. Therefore, combination of products has been rated as «MR unsafe».

Under the examination conditions defined above, it can be assumed that a Ligamys implant leads to a maximum temperature increase of less than 4.4 °C after 15 minutes of continuous scanning. In non-clinical testing, the image artifact caused by the implant extends to about 55mm from the Ligamys, if imaging is carried out using a gradient echo sequence and a 3.0 T MRI system. The magnetically induced displacement force was examined in non-clinical tests. At 3.0 T, a maximum deflection angle of less than 68° was found.

The magnetically induced displacement torque was examined in non-clinical tests. No torque was measured.

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Manufacturer: Mathys Ltd Bettlach